


MK-3CN HEMATOLOGY CONTROL+CRP

2.5mL


ASSAY SHEET

LOT	MK3C232N
	2023-05-05
OPEN	14DAYS


MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	7.67	± 0.80
RBC	10 ⁶ /μL	4.58	± 0.20
HGB	g/dL	13.69	± 0.50
HCT	%	45.5	± 3.0
MCV	fL	99.4	± 7.0
MCH	pg	29.9	± 2.4
MCHC	g/dL	30.1	± 3.0
PLT	10 ³ /μL	308.4	± 50.0
MPV	fL	9.0	± 2.0
RDW-CV	%	15.2	± 3.0
LY%	%	34.20	± 10.00
MO%	%	2.83	± 5.00
GR%	%	62.97	± 10.00
LY	10 ³ /μL	2.62	± 1.00
MO	10 ³ /μL	0.22	± 0.40
GR	10 ³ /μL	4.83	± 1.00
CRP	Precision	≤ 10.0% (20.40-38.00 mg/L)	

Manual Registration

Lot	232
Level	Normal
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU

TOKYO 161-8560 JAPAN

IVD
CE
 0123

0694-901766A