

Bluetooth® Module

QI-832V

General

The QI-832V *Bluetooth*® module is for the following units.

- TEC-5601, TEC-5611, TEC-5621 and TEC-5631 defibrillators
- TEC-8321K, TEC-8322K, TEC-8332K, TEC-8342K and TEC-8352K defibrillators

It lets you send 12 lead analysis results (for TEC-8300 series defibrillator only), waveform reports and self test results to a device such as locally purchased personal computer.

Safety Information

WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

WARNING

Connect the *Bluetooth*® module to the specified defibrillator and follow the specified procedure. Failure to follow this warning may result in electrical shock or injury to the patient and operator.

CAUTION

Before installing or removing the module, make sure that the defibrillator is turned off and the AC/DC module and the battery pack are removed. Otherwise, the operator may receive electrical shock or injury.

CAUTION

Before installing or removing the module, make sure that the defibrillator is turned off and the battery pack is removed. Otherwise, the operator may receive electrical shock or injury.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.



Manufacturer

NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>



European Representative

NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599